

# 2011 Orthodontic Technique Lab Course


## Appliance Projects



**Study Model**



**Lower Invisible Retainer with Reset**



**Upper Wraparound & 'S/Z' Spring**



**Upper Molar Distalizer**



**Distal Jet**



**Pendulum**

Upper 6s = RMO-21  
Upper 4s = RMO-K+



**Upper Schwarz**




**Lower Schwarz**



Upper 6s = RMO-22

**Upper Quadhelix Expander**



**Upper Retainer with Pontic**



**Lower 3x3 Bonded Indirect Lingual Arch**

# 2011 Orthodontic Technique Lab Course




**Indirect Bonding - Dual Laminate Biostar Trays**



Upper 6s = RMO-25

**Upper Transpalatal Arch [TPA]**



Lower Lingual Arch  
Mand 6s = RMO-13

**Lower Lingual 6x6 Arch**



**Haas RME**



Upper 6s = RMO-22  
Upper 4s = RMO-L+



**Upper / Lower Cantilever Herbst**



ORMCO Size - 5  
1<sup>st</sup> Molars [4]



**Finger Spring with Soldered Bow**



**Upper ACCO/Cetlin Molar Uprighter**



**Upper Mouthguard Over Brackets**




# 2011 Orthodontic Technique Lab Course



**Snowman Exercise**



**Scallop Wire Exercise**



**Wirebending & Soldering Exercises**