

Are you considering using a series of expensive aligners for anterior correction?



A Great Lakes Exclusive!

The **NEW** Invisible Inman is a simpler solution!

Moves teeth farther, faster.

The Invisible Inman uses NiTi coil springs to achieve up to 3mm of continuous and controlled lingual-to-labial movement. The piston-like forces of the lingual component gently and precisely move teeth. A labial acrylic shell works to control incisor alignment as teeth are pushed outward. **In the majority of cases, alignment is accomplished using one appliance.**

Invisible for patients. Comfortable, too.

Esthetically pleasing, the appliance has a clear labial component and a lingual wire framework. A low, thin profile makes the Invisible Inman more comfortable than most spring aligners.

Minimal adjustment. Maximum retention.

The NiTi coil spring mechanism retracts to allow you to easily insert and seat the appliance. Adjustments are virtually unnecessary. A thin, flexible acrylic covers the posterior dentition to enhance retention.



Great Lakes
ORTHODONTICS, LTD.
An Employee Owned Company

Photos courtesy of Vincent Bingo, DDS, Honeoye, NY



Patient model before treatment



Model with appliance seated



Progress after 5 months

SM-F-04-10-06

1-800-828-7626 • Fax: 716-871-0550 • www.greatlakesortho.com