



Great Lakes Herbst Advancement Kit

265-107

Stainless steel tubing segments in 1 to 5mm lengths used to advance the Herbst appliance when the desired position can no longer be achieved. Includes: 10 pieces 1mm, 2mm, 3mm, 4mm, and 5mm tubing.

Replacement Tubing:

1mm (10/pkg.)	090-068
2mm (10/pkg.)	090-069
3mm (10/pkg.)	090-070
4mm (10/pkg.)	090-071
5mm (10/pkg.)	090-072

IST-Appliance Kit (Scheu) (Intraoral Snoring Appliance)


Fabricate an easily adjustable snoring appliance using a pressure molding machine. Telescopic hardware controls protrusion, permitting gradual advancement in .25mm increments up to 8mm. Includes: telescopes, acrylic disks, acrylic, fasteners, activation key, retainer box, and instructions.



Kit	Min. Length	Max. Length	Stroke	Activation Length	
HF-9	22mm	39mm	9mm	8mm	191-032
HF-15	28mm	51mm	15mm	8mm	191-033

Telescopic Herbst Hardware (Scheu)

This unique hardware allows the clinician to advance the mandible in .25mm increments (up to 8mm) by making one full turn of the protrusion nut incorporated into the hardware. Also, an end stop in the guiding sleeve prevents the rod and tube from disengaging, reducing the problems associated with loose intraoral hardware.




Important: Herbst Telescopic Hardware should be used for snoring and sleep apnea appliances only.

Components

Short Arm (9mm Arm Length)	090-074
Long Arm (15mm Arm Length)	090-075
Heavy Duty Herbst Sleep Hardware (No Stopper, Telescopic, 2 arms)	090-058
Hex Head Screws and Pivots (4/pkg.) with Allen Wrench	090-077
Attachments with O-rings (4/pkg.)	090-078
IST Positioning Support Arm	090-086
O-rings for IST Appliance (10/pkg.)	090-087

Fixing Elements with Allen Screws and Wrench (4/pkg.)	090-088
Allen Screw for IST Appliance (4/pkg.)	090-089
Hinge System with Retention Arm (2 Herbst Hinges, 4 Attachments, 4 Slotted Screws)	090-090
Slotted Screws (4/pkg.)	090-048
Advancing Keys (25/pkg.)	140-001

Herbst Hardware (Scheu)



Hex Head Kit 090-083

Hex head screw increases accessibility making intraoral manipulation easier and safer. Positive allen wrench insertion prevents screw head from burring for increased patient comfort.


Includes: 2 tubes, 2 plungers, 4 screws, 4 pivots, and screwdriver.

Slot Head Kit 090-050

Includes: 2 tubes, 2 plungers, 4 screws, and 4 pivots. (screwdriver not included)

Pivots & Slot Head Screws (4/pkg.) **090-048**
(090-050 and 090-083 screws are interchangeable)

Herbst Hardware (Dentaurum)



Hex Head Kit 090-045

Hex head screw increases accessibility making intraoral manipulation easier and safer. Positive allen wrench insertion prevents screw head from burring for increased patient comfort.

Includes: 2 tubes, 2 plungers, 4 screws, 4 pivots, and screwdriver.

Slot Head Kit 090-046

Includes: 2 tubes, 2 plungers, 4 screws, and 4 pivots. (screwdriver not included)

Pivots & Slot Head Screws (each) **090-047**

Herbst Components (for standard hardware)

Herbst Plunger (2/pkg.)	090-079
Right Tube	090-080
Left Tube	090-081
Hex Head Screw	090-082
Screwdriver (for slot-head screws)	090-049