

Takes the guesswork out of oral appliance therapy

NEW!





Apnea Guard® Combo Kit

Predicts your patient's therapeutic position

- Precursor to Herbst™* Sleep & other appliances
- Ideal for immediate placement
- Settings can be precisely transferred to a custom device
- Advances at **1mm** increments across full range of protrusion settings
- Uses PVS impression material common to most dental offices
- Produces outcomes equivalent to a custom device

Advances mandible for up to **29** nights

Accommodates full range of arch widths with **3** vertical sizes



255-065

Combo Kit includes:

1-Low (5.5mm), 3-Medium (6.5mm), and 1-High (8mm) Apnea Guards

Replacement Part:

Apnea Guard® Low 255-066 (ea.)



Great Lakes
ORTHODONTICS, LTD.
An Employee Owned Company

*Herbst is a registered trademark of Dentaaurum, Inc.