

**Accurate models are produced every time**

# AquaSpense™



*more reasons. one source*



**Water only dispenser**

### **Specifications:**

Dimensions: 11.75"(W) x 23.25"(H) x 13.12"(D)  
Voltage: 110-240 VAC 50/60Hz  
Amperage: 1.0A  
Weight: 20.5 lbs.

- **Designed primarily for model production**
- **Ensures proper water/powder ratios**
- **Automatically calculates powder weight**
- **Dispenses the correct amount of distilled water**

**AquaSpense takes the human error out of the process to achieve more predictable hardness, expansion working times, and set times.**

**800.828.7626 • Fax: 716.871.0550 • [www.greatlakesortho.com](http://www.greatlakesortho.com)**