



**NEW!**

# Extra Anterior Protection

## Bioplast<sup>®</sup> Xtreme Mouthguard Material

Integrates two degrees of hardness in one material!




- *Transparent, rigid, 'V-shaped' section provides extra anterior protection*
- *You control the number of teeth you want reinforced by changing the position of the model on the 'V'*
- *Increased safety and absorption during hard and flat impacts*

Available colors and sizes:  
10/pkg.



**Xtreme Red:**  
030-090 4mm x 125mm round  
030-093 5mm x 125mm round



**Xtreme Black:**  
030-092 4mm x 125mm round  
030-095 5mm x 125mm round



**Xtreme Blue:**  
030-091 4mm x 125mm round  
030-094 5mm x 125mm round



**Great Lakes**  
ORTHODONTICS, LTD.  
*An Employee Owned Company*

SMBP-339Rev090611  
S-378

**Fabricate your own Xtreme mouthguard, or let Great Lakes Laboratory do it for you**



**NEW!**

## Extra Anterior Protection


# Bioplast<sup>®</sup> Xtreme Mouthguards

**Integrates two degrees of hardness in one material!**

Now you can provide extra anterior protection for your sports patients. Our mouthguards can also be fabricated using Bioplast Xtreme, an innovative new thermal-forming material that integrates two degrees of hardness in a one-step process, single thickness disc.

A transparent, rigid, 'V-shaped' section provides extra anterior protection and absorbs hard and flat impacts. By varying the position of the model on this 'V', we can adjust the number of teeth you want reinforced.

Standard 'V' coverage: Anterior teeth, slightly on the facial of canines.



### Available colors and thicknesses:

- Xtreme Red
- Xtreme Blue
- Xtreme Black
- Medium Weight (4mm)
- Heavy Weight (5mm)

### Additional Options:

- Can be worn over braces
- Breakaway headstrap
- Occlusal indexing

Note: Due to the one-step process, decals and names cannot be added



### For Technical Assistance Contact:

Great Lakes Laboratory Customer Service  
Monday - Friday: 8:00am-5:00pm (EST)

800.828.7626 (US & Canada)

716.871.1161 (Worldwide)

Email: [labtech@greatlakesortho.com](mailto:labtech@greatlakesortho.com)

[www.greatlakesortho.com](http://www.greatlakesortho.com)



**Great Lakes**

ORTHODONTICS, LTD.  
An Employee Owned Company

**Great Lakes Laboratory will fabricate your mouthguard, or you can do it yourself**



800.828.7626 • Fax: 716.871.0550 • [www.greatlakesortho.com](http://www.greatlakesortho.com)