

Innovative. Precise. Portable.

Odyssey[®] Navigator[™] 3Watt Diode Laser



Choose from 30 pre-set clinical procedures

Touch screen monitor inverts for readability in any direction

Rechargeable battery and docking cradle

290-050

Easy snap on uni-dose fiber tips:
 . sterilized
 . pre-cleaved
 . pre-stripped

Ergonomic handpiece design with removable autoclavable sleeve

Benefits for the Doctor:

- Precise tissue modification
- Early bracket placement
- Fast, effective hemostasis
- Minimal or no bleeding
- Reduced treatment time
- Improved productivity
- High patient acceptance and satisfaction
- Increased revenue
- Ability to expand current services

Benefits for the Patient:

- Maximum comfort
- Rapid healing time
- Minimal or no anesthesia
- Less chair time
- Optimal outcome

Additional Features:

- 810 ±20 nm wavelength
- Wireless foot pedal
- Continuous or pulse wave
- 'Quick Set' option
- Unit Size: 7" x 4" x 6"
- Weight: 2.5 lbs.

Applications

Orthodontic

- Crown exposure to aid eruption
- Operculectomy for faster band placement
- Treatment of gingival hyperplasia
- Gingival sculpting
- Ulcer treatment
- Frenectomy

Palliative

- Fibroma removal
- Aphthous ulcer and herpetic lesion treatment
- Bacterial reduction in periodontal pockets

Restorative/Esthetic

- Gingival troughing for taking impressions (no cord to pack)
- Gingival sculpting
- Implant exposure
- Gingivectomy



Before treatment



Partially erupted canine

After treatment



Bracket bonded after tissue sculpting



Gingival hyperplasia



Healthy gingival contour achieved

SMPP99Rev020408

Odyssey® Navigator™ 3Watt Diode Laser



Included with your Odyssey Navigator 3Watt Laser:

- 290-004 Protective glasses - Solid sides style [Medium]
- 290-005 Protective glasses - Plain sides style [Small]
- 290-006 Protective glasses - Window sides style [Large]
- 290-061 Handpiece v2 - 4' length*
- 290-063 Handpiece v2 - Sleeve*
- 290-052 Pack of [50] Standard length (6mm) 400 micron fiber tips
- 290-056 Docking cradle
- 290-057 Power supply
- 290-058 Wireless foot pedal
- 290-059 Lithium polymer battery
- Owner's Manual, Procedure Guide, and Warranty Card

* (2) Handpiece and sleeve assemblies included with initial purchase

Additional accessories:

- 290-051 Pack of [50] Non-contact length (3.5mm) 400 micron fiber tips
- 290-060 Pack of [50] Perio length (9mm) 400 micron fiber tips
- 290-062 Handpiece v2 - 6' length



Great Lakes

ORTHODONTICS, LTD.

An Employee Owned Company