


NEW!

Optimum Reflectivity

doctoreyes™ Ultrabright Dental Mirror

The best possible combination of reflection surface, safe handling, and comfort



191-040 Occlusal Mirror #13 (70mm)

Vaporized glass surface coating is resistant to most acids, alkalines, and chemicals

Advantages:

- *Highly defined intraoral images*
- *True color representation*
- *Great durability for repetitive use*
- *Ultrabright surface coating prevents ghost image reflections*
- *Removable silicone grip handles*
- *Rounded edges for patient comfort*
- *Grips & glass are autoclavable up to 390°F*



Available in lingual, buccal, and various occlusal sizes



Great Lakes
ORTHODONTICS, LTD.
An Employee Owned Company