

Great Lakes

PHILIPS

HEARTSTART

OnSite Defibrillator

HEARTSTART -
designed for an
ordinary person in
an extraordinary
situation

According to the American Heart Association, sudden cardiac arrest (SCA) claims the lives of more than 340,000 people each year in the U.S. alone. Today, fewer than five percent of SCA victims survive, largely because defibrillators do not reach them in time.

Scenario: Someone in your office suddenly collapses and remains unresponsive. Breathing and pulse are absent...

If defibrillation is performed up to one minute after collapse, there is a 90% recovery rate; at three minutes after collapse, there is a 70% recovery rate; at seven minutes after collapse, the estimated save rate drops to 20%. If emergency medical services arrive ten minutes after collapse, the recovery rate is negligible.

- Defibrillators provide the only effective treatment for sudden cardiac arrest, a leading cause of death in the United States.
- Philips, the leader in defibrillation technology, designed the HeartStart OnSite Defibrillator so that a *minimally trained* lay responder can help save a life.
- Defibrillators are fast becoming a mandatory requirement in public buildings, office buildings, schools, and community centers.
- In dental offices, inclusion of a defibrillator is becoming standard and in some states, mandated by law.



Great Lakes
ORTHODONTICS, LTD.
An Employee Owned Company

1-800-828-7626
Fax: 716-871-0550
www.greatlakesortho.com



\$1504.⁰⁰

Includes defibrillator, a set of adult pads, battery, Quick Reference, & instructions.


- Easy to Use
- Lightweight & Mobile
- Safe for Adults & Children

HEARTSTART

OnSite Defibrillator

Compact and Easy to Use. Weighing only 3.3 lbs., the HeartStart defibrillator can be carried to the patient's side. Using clear, natural voice instructions, the OnSite Defibrillator guides the responder through each step of defibrillation and CPR, even providing a reminder to call for emergency help.

Integrated SMART Pad technology senses a responder's actions and automatically adjusts the pace of the voice instructions to match the responder's needs. The HeartStart OnSite Defibrillator can be used on infants and children with specially designed infant/child pads that reduce the amount of energy to an appropriate level.



The American Heart Association estimates that an additional 40,000 lives could be saved each year with improved access to defibrillators.

Ordering the HeartStart Unit and Accessories is easy...

237-001	Philips HeartStart OnSite Defibrillator	\$1504.00
237-002	Battery for HeartStart Defibrillator	153.00
237-003	Adult SMART Pads Cartridge	61.00
237-004	Infant/Child SMART Pads Cartridge	98.00
237-005	Standard Carry Case for HeartStart	134.00
237-006	Slim Carry Case for HeartStart	113.00
237-007	Defibrillator Cabinet, Basic	236.00
237-008	Defibrillator Cabinet, Semi-Recessed	438.00
237-009	Defibrillator Cabinet, Wall Surface	411.00
237-010	Wall Mount Bracket	92.00
237-011	Fast Response Kit	42.00
237-012	Quick Reference, HeartStart	5.00

"Prices subject to change without notice"

1-800-828-7626

Fax: 716-871-0550

www.greatlakesortho.com

Call today for more information or to order the Philips HeartStart OnSite Defibrillator

