

**NEW! From Great Lakes**

**1.3mm**

# Invisacryl **Hard/Soft**

Invisible Retainer Material

***Excellent adaptation  
for an  
ultra comfortable fit***

- *Dual laminate PETG and EVA material*
- *Eliminates drying time prior to thermal forming*
- *Comfortable for patients*
- *Saves time for labs*



**021-061**

10/pkg.  
1.3mm x 125mm round



800.828.7626 • Fax: 716.871.0550 • [www.greatlakesortho.com](http://www.greatlakesortho.com)