

Compact, quiet and versatile

Osada Electric Handpiece



- *Minimal Vibration*
- *Dual Rotation*
- *Keyless Chuck*
- *Easy to Clean*
- *Features 'Instant Stop'*

150-035

- *This handpiece is ideal for lab or chairside use*
- *Designed for both right- and left-handed technicians*
- *Operates with variable speed foot pedal or manual activation*

Specifications

Box: Weight & Dimensions	5 1/4 lb. 4" W x 7" D x 4 1/4" H
Frequency	50/60 Hz
Current Consumption	0.25A & 40W
Variable Rotational Speed	900 - 30,000 rpm
Torque	4.8 N-cm
Variable Micromotor Voltage	DC 2.0 - 48 Volts
Rubber Handpiece Holder	36 x 155 x 16mm 40g
Removable Chucks	Standard 2.35 mm, 1.6 mm or 3.0 mm
Bur Shank Sizes & Length	Length 40 - 50mm

