

Versatile. Easy to use. Compact.

Great Lakes Tack Welder

Versatile

- Ideal for temporary positioning of components prior to soldering or spot welding
- Six power settings for welding a wide range of materials

Easy to use

- Foot pedal permits "hands-free" operation
- Remote cables plug into front of unit for more convenient welding directly on model

Compact

- Vertical profile requires minimal bench space
- Dimensions: 5 1/2" (H) x 2 1/4" (W) x 9" (L)



Efficient design • Requires minimal bench space

240 030 Tack Leads
240-031 Carbon Tip for Auxilliary
Cables (2/pkg.)



Great Lakes

ORTHODONTICS, LTD.
An Employee Owned Company

800.828.7626 • Fax: 716.871.0550 • www.greatlakesortho.com