


**For cutting, trimming, and polishing**

# **Thermal Material Finishing Kit**



085-055

*The same finishing kit that Great Lakes Laboratory uses*

### **Kit includes 8 durable Burs:**

• *Spiral • Triangular • Fine Triangular • Square*  
*For cutting all hard pressure-molded materials*

• *Blue • Black*  
*For trimming and polishing soft materials*

• *White*  
*For final polishing of soft materials*

• *Dimo Pro Finishing Wheel*  
*For final trimming and polishing of hard and soft materials*

