

# Variable temperature control

## Whip Mix Water Bath



190-006 Water Bath  
190-007 Optional Lift Tray

*This thermostatically controlled unit maintains temperatures in a range of 90°F - 160°F to soften a variety of wax and compound materials.*

### Specifications:

Voltage: 115 VAC/60Hz  
Amperage: 1.5A  
Wattage: 175W  
Dimensions: 7.5"(W) x 5.0"(H) x 8.75"(D)  
Temperature: 90°F - 160°F  
Weight: 5 lbs.  
Pan Capacity: 6.5"(W) x 4.25"(H) x .75"(D)

### Features

- Easy to clean exterior
- Removable stainless steel pan for easy sterilization
- Non-mercury thermometer for accurate temperature verification
- Unit stops heating when pan is removed
- Compatible with Hanau Water Bath pans
- UL and CUL approved



**Great Lakes**

ORTHODONTICS, LTD.  
An Employee Owned Company