



PLEASE PRINT

Account # LO _____ **PO #** _____

B PRACTICE TYPE: _____
(i.e., ortho, GP, pedo, prosthodont, oral surgeon, commer. lab)

I DOCTOR: _____

L ADDRESS: _____
(Specify if ship to address is different)

I CITY: _____ STATE: _____ ZIP: _____

N PHONE: (____) _____

G *Provide Country & City Code* FAX: (____) _____

EMAIL: _____

PATIENT: _____ AGE: _____

PLEASE PRINT

DATE SHIPPED: _____

DATE DUE: _____

1 day before appointment

LAB USE ONLY Incoming # cases 1 2 3 4 5+

Customer Used: GLO Acct 2 Day On Call

Portal Upload - No Frt (99) Cust Acct - No Frt (99)

Disinfected 0 1 2 3 4 5 6 7 8 9

Rcvd: _____

B# _____ Via: _____

QC: _____ LPD/Shp: _____

Needs DD Call Rec: _____

NO BITE / MDL - B / C

Doc# _____ Digital ID# _____

Please Provide: Boxes Labels

Rx: _____ Qty: _____
(specify appliance type)

Design Your Own OR Choose a New! eXceed Kit

eXceed Rx- brackets are placed on the pretreatment models using a standard placement method. Bracket placement can be modified using the eXceed software and approved by the doctor.

eXceed Tx- brackets custom positioned using a simulated final setup and virtual arch wire. The setup is based upon clinical instructions from the customer. Bracket positions are selected to optimize this treatment and can be modified by the customer using the eXceed software. Patient photos and Xrays are required.

- Bracket Information:**
- Brackets Enclosed
 - Bracket Inventory at the Lab
 - Use GLO house brackets
 - EasyClip+
 - BioTru Classic
 - BioTru Ceramic
 - BioTru Sapphire
 - Add Molar Tubes

Slot Size: .018 or .022

Prescription: Roth MBT Damon (EasyClip+ Only)

- EasyClip+ Kit**
eXceed Rx + EasyClip+ brackets + BioTru Contour Nitinol Archwires
- BioTru Classic Kit**
eXceed Rx + BioTru Classic twin brackets + BioTru Nitinol Archwires
- BioTru Ceramic Kit**
eXceed Rx + BioTru Ceramic brackets + BioTru Nitinol Archwires
- Aesthetic Express Kit**
eXceed Rx + BioTru Ceramic brackets + Aesthetic Short Term Archwires

Please send the following supplies:

- Reliance Flow-Tain Adhesive (065-040)
- Reliance Assure Universal Bonding Resin (065-011)
- Reliance Liquid Etching Agent (065-004)
- Elastics _____

Tray Fabrication Options:

Dual Tray (hard outer shell, soft inner shell) for all orders. 1 Piece (Standard) Midline Split 3 Pieces Per Arch

License #: _____

Dr. Signature: _____

Master Rx on File # _____

Special Instructions: _____
