



PLEASE PRINT

Account # LO PO # _____

B PRACTICE TYPE: _____
(i.e., ortho, GP, pedo, prosth, oral surgeon, commer. lab)

I DOCTOR: _____

L ADDRESS: _____
(Specify if ship to address is different)

A CITY: _____ STATE: _____ ZIP: _____

D PHONE: (____) _____

R FAX: (____) _____

S EMAIL: _____

PATIENT: _____ AGE: _____

PLEASE PRINT

DATE SHIPPED: _____

DATE DUE: _____
1 day before appointment

LAB USE ONLY Incoming # cases 1 2 3 4 5+

Customer Used: GLO Acct 2 Day On Call

Portal Upload - No Frt (99) Cust Acct - No Frt (99)

Disinfected 0 1 2 3 4 5 6 7 8 9

Rcvd: _____

B# _____ Via: _____

QC: _____ LPD/Shp: _____

Needs DD Call Rec: _____

NO BITE / MDL - B / C

Doc# _____ Digital ID# _____

Please Provide: Boxes Labels

Rx: _____ Qty: _____
(specify appliance type)

Appliance Protection Program (additional fee)

IMPORTANT! Always retain models and bite until appliance is seated.

Damage to models may occur during fabrication, please mark Rx if duplication (additional fee) of model(s) is required.

SMART MOVES ALIGNER

Smart Moves can include a maximum of 3 setup stages that include 2 appliances per stage. In sequence, the first aligner is Hard/Soft followed by a second of hard material. Each should be worn for 2-3 weeks. If further correction is required to reach desired alignment a new impression needs to be provided after initial movement has been achieved. **Please note: when requested, resetting canines may not produce predictable results.**

Appliance Options Upper Lower Both (Please specify)

If not noted, the lab will carve brackets and remove lingual retainers when present.

Brackets or Lingual Retainer: Remain Please Carve

Smart Moves Method, lab determines set up as needed

- max of 3 set up's and 6 trays per impression.

No quote needed, ok to fabricate

Call with quote, before fabrication

2 Set ups, 4 trays

1 Set up, 2 trays

To align reset teeth properly, interproximal reduction may be necessary.

Stripping Preference(s):

Strip and notify me where (Standard)

Strip where indicated UPPER- R 3 2 1 1 2 3 L

Do not strip LOWER- 3 2 1 1 2 3 L

Reset Preference(s):

Reset as needed





Reset where indicated UPPER- R 2 1 / 1 2 L

Over correct where indicated LOWER- 2 1 / 1 2 L

Do not reset

Auxiliaries: **Transfer tray(s) will be provided

Add retention buttons as needed



INVISIBLE RETAINER

Please note: when requested, resetting canines may not produce predictable results.

Appliance Options Upper Lower Both (Please specify)

If not noted, the lab will carve brackets and remove lingual retainers when present.

Brackets or Lingual Retainer: Remain Please Carve

Essix™ Design .75mm (3 x 3 coverage only)

Invisible Retainer 1mm (Full arch coverage)

Scallop Trim Marginal Trim

To align reset teeth properly, it may be necessary to strip interproximally, in these cases:

Stripping Preference(s):

Strip and notify me where

Strip where indicated UPPER- R 3 2 1 1 2 3 L

Do not strip LOWER- 3 2 1 1 2 3 L

Reset Preference(s):

Reset as needed



Reset where indicated UPPER- R 2 1 / 1 2 L

Do not reset LOWER- 2 1 / 1 2 L

Auxiliaries: Plastic Pontic;

Manufacturer Name _____

Shade # _____



Lab Use Only: Block Out Base Dup D.C. Art. Pour
 Resets Pontics Sold. St. M.G. Med. M.G. H.M.G.
Bleaching Tray- Reg. Dura-Soft (Lemon)

License #: _____

Dr. Signature: _____

Master Rx on File # _____

Special Instructions: _____